

2073-A: MORPHOMETRIC DIGITAL PATHOLOGY ANALYSIS REVEALS DIFFERENTIAL EFFECTS OF LANIFIBRANOR AND SEMAGLUTIDE IN THE BIOPSY-CONFIRMED GAN DIO-MASH MOUSE MODEL

Authors:

Cindy Serdjebi¹, Florine Chandès¹, Bastien Lepoivre¹, Susanne E. Pors², Michael Feigh²

¹Biocellvia, Marseille, France
²Gubra A/S, Hørsholm, Denmark

Cindy Serdjebi, R&D Director,
cindy.serdjebi@biocellvia.com

BACKGROUND AND OBJECTIVES

Metabolic dysfunction-associated steatohepatitis (MASH) predisposes to development of advanced fibrosis/cirrhosis. Many clinical trials are ongoing to obtain either significant resolution of MASH without worsening of fibrosis or improvement of fibrosis without worsening of MASH. Semaglutide (glucagon-like-receptor (GLP)-1 agonist) and lanifibranor (pan-peroxisome proliferator activated receptor agonist) are currently in late-stage clinical testing. The present study aimed at investigating the effects of these two monotherapies in the GAN diet-induced obese (DIO) and biopsy-confirmed mouse model of advanced fibrosing MASH using morphometric digital pathology.


MATERIALS AND METHODS

GAN DIO-MASH mice were treated with either vehicle, lanifibranor or semaglutide, and lean chow-fed animals served as control group. The pre- and post-treatment liver biopsies were stained with picosirius red (PSR) and hematoxylin and eosin (H&E) and scanned at the magnification of X20. Histopathological NAFLD Activity Score (NAS) and fibrosis stage were evaluated by Gubra Histopathological Objective Scoring Technique (GHOST) AI-deep learning-based image analysis. In parallel, MorphoQuant, a fully automated and deterministic artificial intelligence assessed steatosis, fibrosis and collagen fiber dimensions (length and width) from PSR-stained sections. Effects of treatments were compared.

CONCLUSION

The GAN DIO-MASH model is highly applicable for profiling novel drug therapies targeting MASH with advanced fibrosis. Notably, morphometric AI-digital pathology showed superior anti-steatotic action for lanifibranor, compared to semaglutide. In addition, the anti-fibrotic effect of lanifibranor, but not semaglutide, was demonstrated, in alignment with the Phase 2 clinical trial data. Importantly, the evaluation of collagen fiber dimensions allows to provide a better understanding of drug effect on fibrosis regression.

.01 STUDY DESIGN AND OUTLINE



Name	Animal Model	Treatment	Subject Number
Chow Vehicle	Chow	Vehicle	9
Vehicle	GAN DIO-MASH	Vehicle	15
Lanifibranor	GAN DIO-MASH	Lanifibranor 30 mg/kg	15
Semaglutide	GAN DIO-MASH	Semaglutide 30 nmol/kg	16

*: $p \leq 0.05$, **: $p \leq 0.01$, ***: $p \leq 0.001$, ****: $p \leq 0.0001$
#: $p \leq 0.05$, ##: $p \leq 0.01$, ###: $p \leq 0.001$, ####: $p \leq 0.0001$

.04 SIGNIFICANT FIBROSIS REGRESSION IN LANIFIBRANOR-TREATED GAN DIO-MASH MICE




Figure 4. Assessment of fibrosis response after 14-weeks treatment with Vehicle, Lanifibranor or Semaglutide in GAN DIO-MASH mice using a Wilcoxon's test A-B-C. Fibrosis score by GHOST for Vehicle-, Lanifibranor- and Semaglutide-treated animals, respectively. D-E-F. Collagen quantification by MorphoQuant for Vehicle-, Lanifibranor- and Semaglutide-treated animals, respectively.

.02 SIGNIFICANT INDUCTION OF STEATOSIS AND FIBROSIS IN GAN DIO-MASH MOUSE MODEL AT BASELINE




Figure 1. Fibrosis and steatosis in pre-treatment liver biopsies from GAN DIO-MASH mice after 50 weeks of diet-induction. A. Fibrosis stage by GHOST. B. Digital quantification of collagen by MorphoQuant. C. Steatosis score by GHOST. D. Steatosis quantification by MorphoQuant. # = t-test between chow vehicle and vehicle, * = Fisher's LSD test between vehicle and treatments

.05 AUTOMATED MORPHOMETRIC ANALYSIS OF FIBROSIS-RELATED ENDPOINTS IN GAN DIO-MASH MICE




Figure 5. Automated morphometric analysis of fibrosis from PSR-stained sections. A. Collagen. B. Collagen width. C. Collagen length..

.03 SIGNIFICANT STEATOSIS REDUCTION IN LANIFIBRANOR- AND SEMAGLUTIDE-TREATED GAN DIO-MASH MICE




Figure 2. Assessment of steatosis response after 14-weeks treatment with Vehicle, Lanifibranor or Semaglutide in GAN DIO-MASH mice using a Wilcoxon's test A-B-C. Steatosis score by GHOST for Vehicle-, Lanifibranor- and Semaglutide-treated animals, respectively. D-E-F. Steatosis quantification by MorphoQuant for Vehicle-, Lanifibranor- and Semaglutide-treated animals, respectively.

.04 AUTOMATED MORPHOMETRIC ANALYSIS OF STEATOSIS-RELATED ENDPOINTS IN GAN DIO-MASH MICE




Figure 3. Automated morphometric analysis of steatosis from PSR-stained sections. A. Steatosis. B. Vesicle area. * = multiple comparison test between vehicle and treatments

.06 MORPHOQUANT ANALYSIS OF COLLAGEN WIDTH IN GAN DIO-MASH MICE




Figure 6. Rainbow color mapped images. A-B. Original image and collagen width recognition for Vehicle-treated animals. C-D. Original image and collagen width recognition for Lanifibranor-treated animals. E-F. Original image and collagen width recognition for Semaglutide-treated animals.

AT A GLANCE

- Lanifibranor and Semaglutide were investigated in GAN DIO-MASH mouse model
- Liver biopsies were assessed for fibrosis and steatosis using GHOST and MorphoQuant (morphometric digital analysis)
- Both Lanifibranor and Semaglutide showed significant decrease of steatosis
- Only Lanifibranor reduced significantly fibrosis, preventing the thickening of collagen fibers
- GAN DIO-MASH mouse model is a highly applicable model for the investigation of NASH treatments with advanced fibrosis
- Evaluation of collagen fiber dimensions allows a better understanding of drug effect